

STEPS TO MAKE A FACE MASK

Protect our health care professionals fighting COVID-19

1 CLEAN YOUR WORK AREA

Using soap and water, followed by a disinfectant, thoroughly clean your work space. Also clean the scissors and your hands. Ensure you don't allow any contaminants into the area and to clean again if you touch your face, phone or any other objects.



2 PREP EQUIPMENT

Set out your roll of Melt-blown Polypropylene, rubber bands and scissors.



3 TEAR ROLL AT PERFERATION

Rip off single sheet of the Melt-blown Polypropylene. This is already folded in half and treated to prevent permiation of the virus.



5 FOLD IN HALF

Fold the entire sheet once more, hot dog style from an inch to a 1/2 inch.



4 FOLD ACCORDION STYLE

Fold back and forth, accordion style, three times, creating four panels that are about one inch each.



7 FOLD OVER FABRIC

Fold the end of the fabric over the rubber band.



6 LOOP RUBBER BAND

Place the rubber band on the end of the fabric about a half inch from the edge.

8 STAPLE

Staple the loose end so that the flat end of the staple is on the inside of the mask. Staple as close to the rubber band as possible



9 REPEAT ON THE OTHER SIDE

Repeat steps 6, 7 and 8 again to complete the otherside of the mask.



10 STORE IN BOX

Put on the first mask and keep on while making the remainder. Place all the completed masks back in the box for storage. Once completed and filled close box and mark with the date and time of completion.

